



IRS, 30.3.21, Labpartner: Dasmih Schleiss Elin Sebesta

- IR spectroscopy has great value in chemistry. Molecules are excited to vibrations and rotations by absorption of photon energy. The vibrations and rotations are specific for certain molecule bonds and their environment. IR spectroscopy can be used for element identification by comparing the obtained spectra to literature references.
- all measurements are made with BRUKER Alpha FT-Spectrometers equipped with ATR and Transmission modules.

→ computer program:
OPUS 7.5



ATR - method



Settings

- Resolution: 2 cm^{-1}
 - Sample scan time: 16 scans (measurement time > 35 s)
 - Background " " " "
 - Save data from 4000 cm^{-1} to: 400 cm^{-1}
 - Result spectrum: Absorbance
 - Spectra extended ATR correction → standard
- ① Background was measured
 - ② Ethanol sample with ATR technique was taken
↳ Cas No.: 64-17-5 → Saved as Ethanol
 - ③ Benzoesäure was measured and then a IR spectrum was taken.
↳ over 99% pure, 1243 g → Saved as Benzoesäure

Test Reproducibility B

- 3 times with Ethanol an IR Spectrum was taken
↳ like ② → LHM from the HCL Shop (not pure)
→ the sample was replaced after the measurements
→ (B1, B2, B3 - Ethanol.dpt)
- 3 times with the Benzoesäure (like from ③)
→ Saved as B1, B2, B3 - Benzoesäure.dpt.

Plastics D → 3 IR Spectra were taken from 3 unknown to said materials

- 1) Handschuh D1 - Handschuh sample was taken
- 2) D2 - ~~sample~~ PET " "
- 3) D3 - ~~CD~~ " "

C DMSO exposed to ambient air

50.3.21, Jasmin
Elin Sebas

- A fresh sample of DMSO was prepared
< 0.03% H_2O
VWR 23500, 7.60

- 10 absorption spectra were measured by
waiting 3 min in between (without changing the sample)
→ saved as C1..10-DMSO.dpt

E → deion. → absolute, > 99.8%

13 H_2O /Ethanol solutions were prepared and an IR
spectrum was taken

- ① 100% H_2O (distilled) 100 μ l H_2O
- ② 100 μ l H_2O 10 μ l Ethanol
- ③ " 20 μ l E
- ④ " 40 μ l E
- ⑤ " 60 μ l E
- ⑥ " 80 μ l E
- ⑦ " 90 μ l E
- ⑧ " 100 μ l E
- ⑨ " 150 μ l E
- ⑩ " 200 μ l E
- ⑪ " 300 μ l E
- ⑫ " ~~300~~ 400 μ l E
- ⑬ 100 μ l Ethanol

→ saved as E1..E13.dpt

→ Samples E-Rum, E-Weisswein, E-Anise
→ 40.0% / → 13.5% / → 35%
→ IR Spectra were taken

50-5-10
Dasmich
Schleiss
Elin ser

F → saved as F-H₂O

1. an IR Spectra from H₂O (deion.) was taken (100 µl)

2. " " D₂O, 99.9% pure (10 µl)
CIL, Inc
DLM4-10